

Алматы (7273)495-231
 Ангарск (3955)60-70-56
 Архангельск (8182)63-90-72
 Астрахань (8512)99-46-04
 Барнаул (3852)73-04-60
 Белгород (4722)40-23-64
 Благовещенск (4162)22-76-07
 Брянск (4832)59-03-52
 Владивосток (423)249-28-31
 Владикавказ (8672)28-90-48
 Владимир (4922) 49-43-18
 Волгоград (844)278-03-48
 Вологда (8172)26-41-59
 Воронеж (473)204-51-73
 Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58
 Иваново (4932)77-34-06
 Иркутск (395)279-98-46
 Казань (843)206-01-48
 Калининград (4012)72-03-81
 Калуга (4842)92-23-67
 Кемерово (3842)65-04-62
 Киров (8332)68-02-04
 Коломна (4966)23-41-49
 Кострома (4942)77-07-48
 Краснодар (861)203-40-90
 Красноярск (391)204-63-61
 Курск (4712)77-13-04
 Курган (3522)50-90-47
 Липецк (4742)52-20-81

Киргизия (996)312-96-26-47

Магнитогорск (3519)55-03-13
 Москва (495)268-04-70
 Мурманск (8152)59-64-93
 Набережные Челны (8552)20-53-41
 Нижний Новгород (831)429-08-12
 Новокузнецк (3843)20-46-81
 Ноябрьск (3496)41-32-12
 Новосибирск (383)227-86-73
 Омск (3812)21-46-40
 Орел (4862)44-53-42
 Оренбург (3532)37-68-04
 Пенза (8412)22-31-16
 Петрозаводск (8142)55-98-37
 Псков (8112)59-10-37
 Пермь (342)205-81-47

Россия (495)268-04-70

Ростов-на-Дону (863)308-18-15
 Рязань (4912)46-61-64
 Самара (846)206-03-16
 Саранск (8342)22-96-24
 Санкт-Петербург (812)309-46-40
 Саратов (845)249-38-78
 Севастополь (8692)22-31-93
 Симферополь (3652)67-13-56
 Смоленск (4812)29-41-54
 Сочи (862)225-72-31
 Ставрополь (8652)20-65-13
 Сыктывкар (8212)25-95-17
 Сургут (3462)77-98-35
 Тамбов (4752)50-40-97
 Тверь (4822)63-31-35

Казахстан (772)734-952-31

Тольяти (8482)63-91-07
 Томск (3822)98-41-53
 Тула (4872)33-79-87
 Тюмень (3452)66-21-18
 Улан-Удэ (3012)59-97-51
 Ульяновск (8422)24-23-59
 Уфа (347)229-48-12
 Хабаровск (4212)92-98-04
 Чебоксары (8352)28-53-07
 Челябинск (351)202-03-61
 Череповец (8202)49-02-64
 Чита (3022)38-34-83
 Якутск (4112)23-90-97
 Ярославль (4852)69-52-93

<https://microlaser.nt-rt.ru> || med@nt-rt.ru

Псевдоколлимированные пучки



Features

- Micron size nearly collimated beams
- Wavelengths from 350nm to 2100nm
- Long depth of focus
- Diffraction limited

Trying to generate a collimated beam that is 100 microns or less in diameter and well behaved is not an easy task. We have simplified this process using our FC3 fiber collimator adjusted to have a nearly collimated beam over some distance.

This collimator is adjusted and set to give a nearly collimated beam that is ten's of microns in diameter at a fixed distance from the collimator exit face. Optically the beam has a long depth of focus, meaning the beam is fairly collimated over a long enough space to be able to insert cuvettes or other samples without worrying about sample position.

Standard distances from the collimator exit face are 100mm or 150mm but other distances can be set as well. We can set to popular wavelength but any wavelength from 350nm to 2100nm can be done as well.

Specifications

FC3 Pseudo Collimated	
Wavelength:	350nm to 2100 nm
Receptacle:	FC or FC/APC
Housing material:	Stainless steel

Ordering Information

Model #	Pseudo collimated Beam size (singlemode fiber)	Working distance (mm)
FC3-405-FC-100	- 90 μ m	100
FC3-450-FC-100		
FC3-488-FC-100		
FC3-520-FC-100		
FC3-635-FC-100		